

FZD9 Antibody



PACO13687

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FZD9

Uniprot

O00144

Synonyms:

frizzled homolog 9 (Drosophila); FZD9; CD349; FZD3 ;

Immunogen:

Human FZD9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A