

MACF1 Antibody



PACO13693

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MACF1

Uniprot

Q9UPN3

Synonyms:

microtubule-actin crosslinking factor 1; MACF1; ABP620; ACF7; FLJ45612; FLJ46776; KIAA0465; KIAA1251; MACF; OFC4 ;

Immunogen:

Human MACF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A