

NPHP4 Antibody



PACO13699

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

Protein Background:

Gene ID:

NPHP4

Uniprot

O75161

Synonyms:

1426aa,158kd;

Immunogen:

Human NPHP4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A