

# OPLAH Antibody



PACO13708

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

OPLAH

**Uniprot**

O14841

**Synonyms:**

5-oxoprolinase (ATP-hydrolysing); OPLAH;5-Opase; DKFZp434H244; OPLA ;

**Immunogen:**

Human OPLAH.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A