

PLXNB2 Antibody



PACO13728

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PLXNB2

Uniprot

O15031

Synonyms:

plexin B2; PLXNB2; KIAA0315; MM1; Nbla00445; PLEXB2; dJ402G11.3 ;

Immunogen:

Human PLXNB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A