

PACO13733

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SCGN

**Uniprot**

O76038

**Synonyms:**

secretagogin, EF-hand calcium binding protein; SCGN; CALBL; DJ501N12.8; SECRET; SEGN; setagin ;

**Immunogen:**

Human SCGN.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A