

THOP1 Antibody



PACO13743

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

THOP1

Uniprot

P52888

Synonyms:

thimet oligopeptidase 1; THOP1; EP24.15; MEPD_HUMAN; MP78; TOP ;

Immunogen:

Human THOP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A