

CD99 Antibody



PACO13757

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CD99

Uniprot

P14209

Synonyms:

CD99 antigen, E2 antigen, Protein MIC2, T-cell surface glycoprotein E2

Immunogen:

Human CD99.

Storage:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Product Images

N/A

N/A