

PACO13917

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:1000-1:10000, IHC:1:50-1:200

Protein Background:

COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The *S. cerevisiae* Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain.

Gene ID:

COX19

Uniprot

Q49B96

Synonyms:

COX19 cytochrome c oxidase assembly factor


Immunogen:

Fusion protein of human COX19.


Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using PACO13917(COX19 Antibody) at dilution 1/50. (Original magnification: x—200).



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using PACO13917(COX19 Antibody) at dilution 1/50. (Original magnification: x—200).