

LGR5 Antibody



PACO14644

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:1000-1:5000, IHC:1:25-1:100

Protein Background:

This gene encodes Orphan receptor. Stem cell marker of the intestinal epithelium and the hair follicle. Target gene of Wnt signaling.

Gene ID:

LGR5

Uniprot

O75473

Synonyms:

leucine-rich repeat containing G protein-coupled receptor 5


Immunogen:

Fusion protein of human LGR5.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using PACO14644(LGR5 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x—200).