

NTMT1 Antibody



PACO16673

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:5000, IHC:1:50-1:200

Protein Background:

The METTL11A gene encodes an N-terminal methyltransferase for the RAN (MIM 601179) guanine nucleotide exchange factor regulator of chromosome condensation 1 (RCC1; MIM 179710). METTL11A enzyme alpha-N-methylates other protein targets such as SET.

Gene ID:

NTMT1

Uniprot

Q9BV86

Synonyms:

N-terminal Xaa-Pro-Lys N-methyltransferase 1


Immunogen:

Fusion protein of human NTMT1.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using PACO16673(NTMT1 Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x—200).