

PACO16947

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:5000, IHC:1:25-1:100

Protein Background:

The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described.

Gene ID:

RASSF4

Uniprot

Q9H2L5

Synonyms:

Ras association (RalGDS/AF-6) domain family member 4


Immunogen:

Fusion protein of human RASSF4.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using PACO16947(RASSF4 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x—200).