

SNX20 Antibody



PACO17140

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:5000, IHC:1:25-1:100

Protein Background:

SNX20 interacts with the cytoplasmic domain of PSGL1 (SELPLG; MIM 600738) and cycles PSGL1 into endosomes.

Gene ID:

SNX20

Uniprot

Q7Z614

Synonyms:

sorting nexin 20


Immunogen:

Fusion protein of human SNX20.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using PACO17140(SNX20 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x—200).