

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:30-1:150

Protein Background:

TRIM 74 (Tripartite motif-containing protein 74) is a possible protein coding regions found at gene location 7q11.23. Tripartite motif (TRIM) proteins play important roles in a variety of cellular functions including cell proliferation, differentiation, development, oncogenesis, and apoptosis. TRIM gene expression analysis in primary human immune cells seem to suggest the involvement of TRIM proteins in also regulating host antiviral activities. The gene encoding TRIM 74 maps to human chromosome 7, which houses over 1,000 genes and comprises nearly 5% of the human genome. Chromosome 7 has been linked to Osteogenesis imperfecta, Pendred syndrome, Lissencephaly, Citrullinemia and Shwachman-Diamond syndrome. The deletion of a portion of the q arm of chromosome 7 is associated with Williams-Beuren syndrome. Deletions of portions of the q arm of chromosome 7 are also seen in a number of myeloid disorders including cases of acute myelogenous leukemia and myelodysplasia.

Gene ID:

TRIM74

Uniprot

Q86UV6

Synonyms:

tripartite motif containing 74


Immunogen:

Full length fusion protein.

Storage:

-20°C; C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using PACO17315(TRIM74 Antibody) at dilution 1/45, on the right is treated with fusion protein. (Original magnification: x—200).