

# FOXD2 Antibody



PACO17888

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:1000-1:2000, IHC:1:15-1:50

**Protein Background:**

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. FOX (Forkhead box) proteins are a family of transcription factors that play important roles in regulating the expression of genes involved in cell growth, proliferation, differentiation, and longevity. Many FOX proteins are important to embryonic development.

**Gene ID:**

FOXD2

**Uniprot**

O60548

**Synonyms:**

forkhead box D2

**Immunogen:**


Synthetic peptide of human FOXD2.

**Storage:**

-20&deg; C, pH7.4 PBS, 0.05% NaN<sub>3</sub>, 40% Glycerol

## Product Images

---



The image on the left is immunohistochemistry of paraffin-embedded Human brain cancer tissue using PACO17888(FOX2 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x—200).