

# TICAM1 Antibody



PACO21086

---

## Product Information

**Size:**

100ul(100ug)

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:50-1:200

**Protein Background:****Gene ID:**

TICAM1

**Uniprot**

Q8IU6

**Synonyms:**

TRIF; PRVTIRB; TICAM-1

**Immunogen:**


Recombinant protein of human TICAM1.

**Storage:**

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of lung cell lysate using TICAM1 antibody.