

GRIA3 Antibody



PACO21090

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

GRIA3

Uniprot

P42263

Synonyms:

GLUR3; GLURC; GluA3; MRX94; GLUR-C; GLUR-K3


Immunogen:

Recombinant protein of human GRIA3.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of brain cell lysate using GRIA3 antibody.