

ADIPOR1 Antibody



PACO21142

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:

Gene ID:

ADIPOR1

Uniprot

Q96A54

Synonyms:

Adiponectin receptor protein 1; Progestin and adipoQ receptor family member 1 ; CGI-45; PAQR1; ADIPOR1


Immunogen:

Recombinant protein of human ADIPOR1.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of Y79 cell line, using ADIPOR1 antibody.