

HDAC9 Antibody



PACO21144

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

HDAC9

Uniprot

Q9UKV0

Synonyms:

HDAC9; DKFZp779K1053; HD7; HD7b; HD9; HDAC; HDAC7; HDAC7B; HDAC9B;
HDAC9FL; HDRP; KIAA0744; MITR


Immunogen:

A Synthetic peptide of human HDAC9.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using HDAC9 antibody.