

SYK Antibody



PACO21256

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

SYK

Uniprot

P43405

Synonyms:

SYK; DKFZp313N1010; FLJ25043; FLJ37489

Immunogen:

Recombinant protein of human SYK.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot analysis of extracts of Raji cell lines, using SYK antibody.

