

ABCC4 Antibody



PACO21292

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:

Gene ID:

ABCC4

Uniprot

O15439

Synonyms:

ABCC4; EST170205; MOAT-B; MOATB; MRP4;

Immunogen:

Recombinant protein of human ABCC4.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A