

IRF3 (Ab-396) Antibody



PACO21334

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

IRF3

Uniprot

Q14653

Synonyms:

interferon regulatory factor 3


Immunogen:

Peptide sequence around aa.394-398(H-I-S-N-S) derived from Human IRF-3.

Storage:

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 and 293 cells using IRF-3 (Ab-396) Antibody.