

# ANXA5 Antibody



PACO21347

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:**

**Gene ID:**

ANXA5

**Uniprot**

P08758

**Synonyms:**

Anchorin CII; Annexin 5; Annexin A5; Annexin V; Annexin-5; Annexin5; AnnexinA5; AnnexinV; ANX 5; ANX A5; ANX5; ANXA5; ANXA5\_HUMAN; Calphobindin I; CBP I; CBP-I; Endonexin II; ENX 2; ENX2; Lipocortin V; PAP I; PAP-I; Placental anticoagulant protein 4;

**Immunogen:**


Purified Recombinant Human ANXA5 protein fragments expressed in E.coli.

**Storage:**

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot detection of ANXA5 antibody in 3T3, HeLa, C6, COS7 and MCF7 cell lysates using ANXA5 antibody (1:2000 diluted). Predicted band size:30kDa. Observed band size:30kDa. Exposure time:20s.