

TP53 (Ab-6) Antibody



PACO21357

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TP53

Uniprot

P04637

Synonyms:

Tumor suppressor p53; Phosphoprotein p53; Antigen NY-CO-13; TP53;


Immunogen:

Peptide sequence around aa.4~8 (P-Q-S-D-P) derived from Human p53.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from MDA cells using p53(Ab-6) Antibody.