

# TP53 (Ab-9) Antibody



PACO21358

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:**

**Gene ID:**

TP53

**Uniprot**

P04637

**Synonyms:**

Tumor suppressor p53; Phosphoprotein p53; Antigen NY-CO-13; TP53;

**Immunogen:**


Peptide sequence around aa.7~11 (D-P-S-V-E) derived from Human p53.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from 293, HT29 and MDA cells using p53(Ab-9) Antibody.