

# TP53 (Ab-18) Antibody



PACO21359

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:****Protein Background:****Gene ID:**

TP53

**Uniprot**

P04637

**Synonyms:**

Tumor suppressor p53; Phosphoprotein p53; Antigen NY-CO-13

**Applications:**

ELISA, WB

**Immunogen:**

Peptide sequence around aa. 16~20 (Q-E-T-F-S) derived from Human p53.

**Recommended dilutions:**


ELISA:1:2000-1:10000, WB:1:500-1:1000

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HT29 cells using p53(Ab-18) Antibody.