

MAPT (Ab-214) Antibody



PACO21369

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

MAPT

Uniprot

P10636

Synonyms:

MAPT; MTAPT; MTBT1; Neurofibrillary tangle protein; PHF-tau

Applications:

ELISA, IF

Immunogen:

Peptide sequence around aa.212~216 (T-P-S-L-P) derived from Human Tau.


Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:200

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunofluorescence staining of methanol-fixed HeLa cells using Tau(Ab-214) Antibody.