

ICAM1 (Ab-512) Antibody




PACO21374

Product Information

Size:	Protein Background:
100ul	
Reactivity:	Gene ID:
Human	ICAM1
Source:	Uniprot
Rabbit	P05362
Isotype:	Synonyms:
IgG	ICA1; ICAM1;
Applications:	Immunogen:
ELISA, WB	Peptide sequence around aa. 510~514 (K-K-Y-R-L) derived from Human ICAM-1.
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:500-1:1000	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from MDA cells using ICAM-1(Ab-512) Antibody.