

RB1 (Ab-795) Antibody



PACO21375

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

RB1

Uniprot

P06400

Synonyms:

P105-RB; PP105; PP110; RB-1; RB1

Applications:

ELISA, WB

Immunogen:

Peptide sequence around aa.793~797 (P-S-S-P-L) derived from Human Rb.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from MOLT cells using Rb(Ab-795) Antibody and the same antibody preincubated with blocking peptide.