

# TYK2 (Ab-1054) Antibody



PACO21378

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:****Gene ID:**

TYK2

**Uniprot**

P29597

**Synonyms:**

tyrosine kinase 2; JTK1;

**Immunogen:**


Peptide sequence around aa.1052~1056 (H-E-Y-Y-R) derived from Human TYK2.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using TYK2(Ab-1054) Antibody(left) or the same antibody preincubated with blocking peptide(right).