

RELA (Ab-505) Antibody



PACO21389

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:**Gene ID:**

RELA

Uniprot

Q04206

Synonyms:

NFKB3; RELA; TF65; Transcription factor p65; p65


Immunogen:

Peptide sequence around aa.503~507 (L-V-T-G-A) derived from Human NFkB-p65.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NFKB-p65(Ab-505) Antibody(left) or the same antibody preincubated with blocking peptide(right).