

MAPK3 (Ab-204) Antibody



PACO21402

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

MAPK3

Uniprot

P27361

Synonyms:

Extracellular signal-regulated kinase 1


Immunogen:

Peptide sequence around aa. 202~206 (T-E-Y-V-A) derived from Human p44/42 MAP Kinase.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from PC12 and 293 cells using p44/42 MAP Kinase(Ab-204) Antibody.