

RELB (Ab-573) Antibody



PACO21404

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

RELB

Uniprot

Q01201

Synonyms:

I-Rel


Immunogen:

Peptide sequence around aa. 550~554 (L-L-S-P-G) derived from Human RelB.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC, 293 and HeLa cells using RelB(Ab-573) Antibody.