

DAB1 (Ab-232) Antibody



PACO21406

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

DAB1

Uniprot

O75553

Synonyms:

Disabled homolog 1

Applications:

ELISA, IHC

Immunogen:

Peptide sequence around aa.230~234 (G-V-Y-D-V) derived from Human Dab1.


Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Dab1(Ab-232) Antibody(left) or the same antibody preincubated with blocking peptide(right).