

# GAB1 (Ab-627) Antibody



PACO21414

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:****Gene ID:**

GAB1

**Uniprot**

Q13480

**Synonyms:**

Growth factor receptor bound protein 2-associated protein 1

**Immunogen:**


Peptide sequence around aa.625~629 (V-E-Y-L-D) derived from Human Gab1.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from A431 and MCF cells using Gab1 (Ab-627) Antibody.