

GSK3A/GSK3B (Ab-279/216) Antibody



PACO21418

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

GSK3A/GSK3B

Uniprot

P49840/P49841

Synonyms:

Factor A; GSK-3 alpha/beta; kinase GSK3-alpha/beta


Immunogen:

Peptide sequence around aa.277~281/214~218 (V-S-Y-I-C) derived from Human GSK3 /.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HL60, HT29 and 293 cells using GSK3- α ; / β (Ab-279/216) Antibody.