

PIK3R1/PIK3R3 (Ab-467/199) Antibody



PACO21428

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

PIK3R1/PIK3R3

Uniprot

P27986/Q92569

Synonyms:

p85, AGM7, p85-ALPHA, p55, p55-GAMMA


Immunogen:

Peptide sequence around aa.465~469(L-Y-E-E-Y) derived from Human PI3 Kinase p85/p55.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using PI3K P85 a/g/b(Ab-467/199/464) Antibody(left) or the same antibody preincubated with blocking peptide(right).