

MTOR (Ab-2481) Antibody



PACO21430

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

MTOR

Uniprot

P42345

Synonyms:

FRAP; FRAP1; FRAP2; RAFT1; Rapamycin target protein

Applications:

ELISA, IHC

Immunogen:

Peptide sequence around aa.2479~2483 (I-H-S-F-I) derived from Human mTOR.


Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using mTOR(Ab-2481) Antibody(left) or the same antibody preincubated with blocking peptide(right).