

# FABP4 Antibody



PACO21477

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000

**Protein Background:**

Lipid transport protein in adipocytes. Binds both long chain fatty acid, and retinoic acid, Delivers long-chain fatty acid, and retinoic acid, to their cognate receptors in the nucleus Marr E. , Tardie M. , Carty M. , Brown Phillips T. , Wang IK. , Soeller W. , Qiu X. , Karam G. Acta Crystallogr. F 62:1058-1060(2006).

**Gene ID:**

FABP4

**Uniprot**

P15090

**Synonyms:**

aP2; ALBP; AFABP; A-FABP;

**Immunogen:**


Peptide sequence around aa.123~127(T-S-T-R-V)derived from Human FABP4.

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extract from 3T3, HeLa cells using FABP4 Antibody.