

GAD2 Antibody



PACO21488

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Catalyzes the production of GABA.

Gene ID:

GAD2

Uniprot

Q05329

Synonyms:

MGC161605; MGC161607;


Immunogen:

Peptide sequence around aa.148~152 (E-L-A-D-Q) derived from Human GAD65 (GAD2).

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from Rat brain, Mouse brain tissue and HeLa, 3T3 cells using GAD65(GAD2) Antibody.