

MELK Antibody



PACO21502

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

MELK

Uniprot

Q14680

Synonyms:

hPK38; KIAA0175; hMELK


Immunogen:

Peptide sequence around aa.626~630 (I-R-R-Q-R) derived from Human MELK.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from Rat brain tissue using MELK Antibody.