

PACO21510

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization. Sugiura T. , Takamatsu H. , Kudo A. , Amann E. Protein Expr. Purif. 6:305-311(1995).

Gene ID:

Postn

Uniprot

Q62009

Synonyms:

POSTN; OSF2; OSF-2; PN; osteoblast specific factor


Immunogen:

Peptide sequence around aa. 125~129(R-E-E-I-E) derived from Mouse Periostin.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from Rat brain and Mouse brain Tissue using Periostin Antibody.