

MAP1LC3A Antibody



PACO21513

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Probably involved in formation of autophagosomal vacuoles (autophagosomes).

Gene ID:

MAP1LC3A

Uniprot

Q9H492

Synonyms:

LC3A; ATG8E; MAP1ALC3; MAP1BLC3;


Immunogen:

Peptide sequence around aa. 112-116(V-Y-A-S-Q) derived from Human LC3.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from K-562, HeLa, JK and 293 cells using LC3 Antibody.