

FOXA2 Antibody



PACO21522

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

FOXA2

Uniprot

Q9Y261

Synonyms:

HNF3B; TCF3B;

Applications:

ELISA, WB

Immunogen:

Peptide sequence around aa.435~439(A-S-P-L-A) derived from Human FoxA2(HNF3b).


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293, HT29, HepG2 cells using FoxA2(HNF3b) Antibody.