

ANXA6 Antibody



PACO21532

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

ANXA6

Uniprot

P08133

Synonyms:

ANX6; CBP68;


Immunogen:

Peptide sequence around aa.6~10(Q-G-A-K-Y) derived from Human Annexin VI.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 and HeLa cells using Annexin VI Antibody.