

# Gabra3 Antibody



PACO21542

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Rat, Mouse

**Source:**

Rabbit

**Isotype:****Protein Background:****Gene ID:**

Gabra3

**Uniprot**

P20236

**Synonyms:**

Gamma-aminobutyric acid, receptor subunit alpha-3; Gabra3;

**Applications:**

ELISA, WB

**Immunogen:**

Peptide sequence around aa. 33~37(R-R-Q-E-P)derived from Rat GABA A Receptor a3.

**Recommended dilutions:**


ELISA:1:2000-1:10000, WB:1:500-1:1000

**Storage:**

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extract from rat brain and mouse brain tissue and C6 cells using GABA A Receptor  $\alpha 3$  Antibody.