

CNP Antibody



PACO21544

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Nucleoside 2',3'-cyclic phosphate + H₂O = nucleoside 2'-phosphate.

Gene ID:

CNP

Uniprot

P09543

Synonyms:

CNP1


Immunogen:

Peptide sequence around aa.80~84 (M-V-S-A-D).

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from Rat brain, Mouse brain and Rat hippocampus Tissue using CNPase Antibody.