

TFDP1 Antibody



PACO21548

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TFDP1

Uniprot

Q14186

Synonyms:

Dp-1; TFDP1; DRTF1


Immunogen:

Peptide sequence around aa.13~17(E-L-K-V-F) derived from Human DP-1.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa, HUVEC and HepG2 cells using DP-1 Antibody.