

E2F2 Antibody



PACO21551

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

E2F2

Uniprot

Q14209

Synonyms:

E2F-2


Immunogen:

Peptide sequence around aa.425~429(L-F-D-S-Y) derived from Human E2F-2.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa and K562 cells using E2F-2 Antibody.