

ENO3 Antibody



PACO21552

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

ENO3

Uniprot

P13929

Synonyms:

MSE; GSD13;


Immunogen:

Peptide sequence around aa.49~53(L-R-D-G-D) derived from Human b-Enolase(ENO-3).

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 3T3, HepG2 and C2C12 cells using β -Enolase(ENO-3) Antibody.